



PRECISION-CUT GUT SLICES

PROTOCOL FOR CUTTING SLICES OF FRESH GUT TISSUE

Key to reading the protocol:

✓ Rationale for procedural step

♠ Tips & Tricks

1. Prepare Krebs-Henseleit Buffer (KHB) for slicing intestine and storing intestinal slices.

✓ Here's the simple recipe for KHB:

KHB (slicing buffer)	Concentration (mM)
NaCl	118
KCl	5
MgSO ₄ ·7H ₂ O	1.1
KH ₂ PO ₄	1.2
NaHCO ₃	25
CaCl ₂ ·2H ₂ O	2.5
D-Glucose	25
HEPES	9

2. Prepare Williams medium E (WME, containing L-glutamine)

✓ Here's the simple recipe for KHB:

WME	Concentration (mM)
WME + L-glutamine	5
D-Glucose	14
Gentamycin	50 µg/ml

Amphotericin B	2.5 $\mu\text{g/ml}$
----------------	----------------------

- √ Amphotericin B is added to inhibit growth of fungi and yeasts.
- 3. Fill the culture plates you will use with the incubation medium.
- 4. Warm and oxygenate culture plates with 95% O₂/5% CO₂ at 37 C for at least 30 min.
- 5. Prepare and keep warm a 2-3% solution of agarose at 37 C.
- 6. To obtain intestinal tissue, anesthetize the animal and retrieve the gut as quickly as possible, ideally under 5 minutes. Remove any fatty tissue and flush the intestines by pipetting ice-cold oxygenated KHB solution. Transfer cleaned pieces of intestinal tissue to ice-cold oxygenated KHB.
- √ Intestinal tissue is very sensitive to ischemia, so the time taken for dissection and incubation of precision-cut gut tissue should be as fast as possible.
- 7. Tighten one end of the gut segment with surgical thread. Using a transfer pipette, fill the intestine with 2-3% agarose. Close the end of intestine with thread or a clamp, then transfer into cold oxygenated KHB to help solidify the agarose into a gel (<1 min required).
- 8. Once agarose has solidified, cut the intestinal segment into sections that you want to slice.
- 9. Glue the tissue sample onto the Compressstome[®] specimen syringe.
- 10. Draw the syringe downward to bring the intestinal tissue core sample into the syringe.
- 11. Fill the syringe with 2% agarose (Sigma A-0701, low gelling point, incubated at ~37°C).
- √ Order a Starter Kit or additional agarose or blades directly from our website at <http://www.precisionary.com/starter-kit!>
- 12. Cool the entire contents of the specimen syringe with the chilling block. The intestinal tissue is now embedded in agarose. The agarose will solidify enough for stable sectioning.
- 13. Load the specimen syringe onto the Compressstome[®] slicer.
- 14. The protocol is complete for preparing the intestinal specimen for sectioning. Proceed from here with normal Compressstome[®] sectioning procedures.
- ♠ What are the optimal settings on the Compressstome[®] for cutting live intestinal slices? Try a speed (Advance) of 4 and an oscillation of 5-7.

References

* Uses the Compressstome[®] for successful eye tissue slices.

1. Buffington SA, Di Prisco GV, Auchtung TA, Ajami NJ, Petrosino JF, Costa-Mattioli M. Microbial Reconstitution Reverses Maternal Diet-Induced Social and Synaptic Deficits in Offspring. *Cell*. 2016 Jun 16;165(7):1762-75.

2. de Graaf IA, Olinga P, de Jager MH, Merema MT, de Kanter R, van de Kerkhof EG, Groothuis GM. Preparation and incubation of precision-cut liver and intestinal slices for application in drug metabolism and toxicity studies. *Nat Protoc.* 2010 Sep;5(9):1540-51.
3. de Kanter R, Tuin A, van de Kerkhof E, Martignoni M, Draaisma AL, de Jager MH, de Graaf IA, Meijer DK, Groothuis GM. A new technique for preparing precision-cut slices from small intestine and colon for drug biotransformation studies. *J Pharmacol Toxicol Methods.* 2005 Jan-Feb;51(1):65-72.
4. de Kanter R, Monshouwer M, Meijer DK, Groothuis GM. Precision-cut organ slices as a tool to study toxicity and metabolism of xenobiotics with special reference to non-hepatic tissues. *Curr Drug Metab.* 2002 Feb;3(1):39-59.