



HEART TISSUE SECTIONING

PROTOCOL FOR CUTTING SLICES OF HEART TISSUE

Key to reading the protocol:

√ Rationale for procedural step

♠ Tips & Tricks

1. Deeply anesthetize animal with pentobarbital, 70 mg/kg.
 - √ Note that you can perfuse the animal in a variety of ways, including gravity-feed, pump, or push-syringe.
2. Remove hearts by dissection and perfuse in Langendorff-mode with bicarbonate-buffered solution.
 - √ Note that fat accumulation at the epicardial surface increases with animal age. Fat tissue is difficult to cut and can blunt the blade, so we recommend careful manual removal of fat tissue before slicing.
3. Load dye via the coronary circulation, apply by injecting into the aortic cannula. Use 22 μ l of a solution containing the voltage-sensitive dye di-4-ANBDQPPQ and then Pluronic F-127 (2 μ l of a 20% stock solution in DMSO. Add over a 4 to 5 minute time period.
4. After dye loading, perfuse hearts at room temperature with BDM-containing HEPES-buffered solution.
5. Select a section of the heart that you would like to take cut for slices.
6. Glue the tissue sample onto the Compressstome® specimen syringe.
7. Draw the syringe downward to bring the brain tissue core sample into the syringe.
8. Fill the syringe with 4% agarose (Sigma A-0701, low gelling point, incubated at \sim 37°C).
 - √ Order a Starter Kit or additional agarose or blades directly from our website at <http://www.precisionary.com/starter-kit> !
9. Cool the entire contents of the specimen syringe with the chilling block. The heart tissue is now embedded in agarose. The agarose will solidify enough for stable sectioning.

10. Load the specimen syringe onto the Compresstome® slicer.
11. The protocol is complete for preparing the heart specimen for sectioning. Proceed from here with normal Compresstome® sectioning procedures.

References

*** Uses the Compresstome® for successful heart tissue slices.**

1. Chang AJ, Ortega FE, Riegler J, Madison DV, Krasnow MA. Oxygen regulation of breathing through an olfactory receptor activated by lactate. *Nature*. 2015 Nov 12;527(7577):240-4. doi: 10.1038/nature15721. PubMed PMID: 26560302; PubMed Central PMCID: PMC4765808.
2. George AA, Lyons-Warren AM, Ma X, Carlson BA. A diversity of synaptic filters are created by temporal summation of excitation and inhibition. *J Neurosci*. 2011 Oct 12;31(41):14721-34. doi: 10.1523/JNEUROSCI.1424-11.2011. PubMed PMID:21994388; PubMed Central PMCID: PMC3214667.
3. Hoffman JJ, Johnson BL, Holland MR, Fedewa RJ, Nair A, Miller JG. Layer-dependent variation in the anisotropy of apparent integrated backscatter from human coronary arteries. *Ultrasound Med Biol*. 2011 Apr;37(4):632-41. doi:10.1016/j.ultrasmedbio.2011.01.007. Epub 2011 Mar 3. PubMed PMID: 21376456;PubMed Central PMCID: PMC3063363.
4. Kamari Y, Cohen H, Shaish A, Bitzur R, Afek A, Shen S, Vainshtein A, Harats D. Characterisation of atherosclerotic lesions with scanning electron microscopy(SEM) of wet tissue. *Diab Vasc Dis Res*. 2008 Mar;5(1):44-7. doi:10.3132/dvdr.2008.008. PubMed PMID: 18398812.
5. Komatsu S, Kitazawa T, Ikebe M. Visualization of stimulus-specific heterogeneous activation of individual vascular smooth muscle cells in aortic tissues. *J Cell Physiol*. 2017 Mar 11. doi: 10.1002/jcp.25903. [Epub ahead of print] PubMed PMID: 28295256.
6. Mukherjee S, Trice J, Shinde P, Willis RE, Pressley TA, Perez-Zoghbi JF. Ca²⁺ oscillations, Ca²⁺ sensitization, and contraction activated by protein kinase C in small airway smooth muscle. *J Gen Physiol*. 2013 Feb;141(2):165-78. doi:10.1085/jgp.201210876. PubMed PMID: 23359281; PubMed Central PMCID: PMC3557311.
7. Wang K, Lee P, Mirams GR, Sarathchandra P, Borg TK, Gavaghan DJ, Kohl P, Bollensdorff C. Cardiac tissue slices: preparation, handling, and successful optical mapping. *Am J Physiol Heart Circ Physiol*. 2015 May 1; 308(9): H1112–H1125.