



PROTOCOL FOR IN-SITU HYBRIDIZATION (DNA & RNA)

Key to reading the protocol:

√ Rationale for procedural step

♠ Tips & Tricks

1. For [in-situ hybridization](#) (ISH), fresh tissue should be fixed in 4% paraformaldehyde overnight, then rinsed in 1X PBS the next day. Tissue should be sectioned with the [Compresstome®](#) for desired slices.
 - √ Be sure to not “over-fix” tissue in paraformaldehyde! Doing so results in excessive cross-linking of proteins that can interfere with getting good signals for in-situ hybridization.
 - ♠ What slice thickness should you aim for? For ISH, slices should be [3 µm to 7 µm thick](#). They can be as thick as 10 µm.
2. Mount free-floating slices onto Superfrost-Plus slides (Fisher Scientific).
 - √ Why use Superfrost-PLUS slides? Because there is a (+) charge on these types of slides and allow tissue slices to adhere to the slide surface. Therefore, when you are processing the slides, the slices won't just fall off.
3. Dry slices on slide using a slide warmer for 30 min at 37°C.
4. Fix slides again in 4% paraformaldehyde in diethyl pyrocarbonate (DEPC) H₂O for 20 min.
5. Rinse slides with DEPC PBS (2 x 5 min).
6. Deproteinase slices with proteinase K treatment (50 µg/ml; 1M Tris HCl at pH 7.5, and 0.1M EDTA) for 30 min at 25°C.
7. After deproteinase, rinse sections in DEPC PBS (2 x 5 min).
8. Fix slides in 4% paraformaldehyde in DEPC H₂O for 20 min.
9. Briefly rinse slides in DEPC H₂O.

10. Next, hybridize slices in hybridization buffer containing the oligonucleotides you are using. Incubate for 24 hours at 37°C in a humid chamber.
11. After hybridization, wash sections in 0.2X SSC (4 X 20 min each) at 47°C.
12. Block slices in blocking protein in Tris buffered saline (TBS) for 30 min (Perkin Elmer Amplification Kit).
13. Incubate slides in digitoxigenin-horseradish peroxidase in TBS containing Tween 20 (TBST; 8 µl/ml) for 30 min at 25°C.
14. Rinse slides in TBST (3 X 5 min).
15. Incubate slides in tyramide solution (1:50 in kit diluent) for 10 min at 25°C.
16. Rinse slides in TBST (3 X 5 min).
17. Incubate slides in streptavidin-F (1:50) in TNB buffer (0.05 g/ml block protein in Tris-saline) for 30 min.
18. Rinse slides in TBST (3 X 5 min).
19. Air dry slides and mount coverslips using Vectashield Mounting Media.

References

*** Uses the Compresstome® for successful in-situ hybridization.**

1. **Abdelaal HM, Kim HO, Wagstaff R, Sawahata R, Southern PJ, Skinner PJ. Comparison of Vibratome and Compresstome sectioning of fresh primate lymphoid and genital tissues for in situ MHC-tetramer and immunofluorescence staining. Biol Proced Online. 2015 Jan 7;17(1):2.**
2. **Duncan J, Kersigo J, Gray B, Fritzscht B. Combining lipophilic dye, in situ hybridization, immunohistochemistry, and histology. J Vis Exp. 2011 Mar 17;(49).**
3. **Li S, Folkvord JM, Rakasz EG, Abdelaal HM, Wagstaff RK, Kovacs KJ, Kim HO, Sawahata R, MaWhinney S, Masopust D, Connick E, Skinner PJ. Simian Immunodeficiency Virus-Producing Cells in Follicles Are Partially Suppressed by CD8+ Cells In Vivo. J Virol. 2016 Nov 28;90(24):11168-11180.**
4. **Lu Y, Ye JH. Glycine-activated chloride currents of neurons freshly isolated from the prefrontal cortex of young rats. Brain Res. 2011 Jun 1;1393:17-22.**
5. **Peça J, Feliciano C, Ting JT, Wang W, Wells MF, Venkatraman TN, Lascola CD, Fu Z, Feng G. Shank3 mutant mice display autistic-like behaviours and striatal dysfunction. Nature. 2011 Apr 28;472(7344):437-42.**
6. **Selever J, Kong JQ, Arenkiel BR. A rapid approach to high-resolution fluorescence imaging in semi-thick brain slices. J Vis Exp. 2011 Jul 26;(53).**

7. Vega-Zuniga T, Marín G, González-Cabrera C, Planitscher E, Hartmann A, Marks V, Mpodozis J, Luksch H. Microconnectomics of the pretectum and ventral thalamus in the chicken (*Gallus gallus*). *J Comp Neurol*. 2016 Aug 1;524(11):2208-29.