



## VF-700-0Z Technical specs

- **Advance Speed:** 0-20 mm / s, adjustable
- **Return Speed:** 20 mm / s, fixed
- **Vibration Frequency:** 0-20 Hz, adjustable
- **Vibration Amplitude:** 0.3 mm, fixed
- **Z-axis Vibration:** ~ 0 $\mu$ m
- **Blade:** double edge razor/ceramic blades
- **Cutting Angle:** 20 degrees fixed (standard)
- **Thickness Adjustment:** digital/automatic
- **Micrometer Resolution:** 1  $\mu$ m
- **Maximum Tissue Diameter:** 22 mm
- **Slice Thickness:**
  - Minimum: 4  $\mu$ m
  - Maximum: 250  $\mu$ m
- **Maximum Tissue Length:** 25 mm
- **Cutting Mode:** single/continuous selectable
- **Bath:** 210 x 80 x 35 mm
- **Power Source:** DC 13-15 V
- **Power Consumption:** 4 W
- **Dimension ( L x W x H ): 300 x 280 x 280 mm**
- **Weight:** 9kg

