

# VF-700-0Z

## Technical specs



**PRECISIONARY**



Advance Speed	0-20 mm / s, adjustable
Return Speed	20 mm / s, fixed
Vibration Frequency	0-20 Hz, adjustable
Vibration Amplitude	0.3 mm, fixed
Z-axis Vibration	~0 $\mu$ m
Blade	double edge razor blade
Cutting Angle	20 degrees, fixed (standard)
Thickness Adjustment	digital adjustment, automatic
Minimum Slice Thickness	10 $\mu$ m
Maximum Adjustable Thickness	250 $\mu$ m
Slice Thickness Resolution	1 $\mu$ m
Maximum Tissue Diameter	20 mm
Maximum Tissue Length	25 mm
Cutting Mode	single/continue selectable
Bath	210 x 80 x 35 mm
Power Source	DC 13-15 V
Power Consumption	4 W
Dimension ( L x W x H )	300 x 280 x 280 mm
Weight	9kg

Precisionary Instruments LLC  
90 Canal Street, Suite 400  
Boston, MA 02114

Tel: 617-682-0586  
Fax: 1(866)424-2217  
Email: [info@precisionary.com](mailto:info@precisionary.com)