


VF-900-0Z

Technical specs



PRECISIONARY



Advance Speed	0-20 mm / s, adjustable
Return Speed	20 mm / s, fixed
Vibration Frequency	0-20 Hz, adjustable
Vibration Amplitude	0.4 mm, fixed
Z-axis Vibration	< 5 μ m
Blade	stainless steel
Cutting Angle	20 degrees, fixed (standard)
Thickness Adjustment	digital/automatic
Slices Thickness	Minimum: 40 μ m Maximum: 500 μ m
Slice Thickness Resolution	2 μ m
Tissue Diameters	170 mm
Maximum Tissue Length	120 mm
Cutting Mode	single/continuous selectable
Bath	730 x 250 x 45 mm
Power Source	DC 13 to 15V
Power Consumption	8W
Dimension	760x 430 x 450 mm
Weight	30 kg

Precisionary Instruments LLC
90 Canal Street, Suite 400
Boston, MA 02114

Tel: 617-682-0586
Fax: 1(866)424-2217
Email: info@precisionary.com