

VF-700-0Z

Technical specs



PRECISIONARY



| | |
|------------------------------|-------------------------------|
| Advance Speed | 0-20 mm / s, adjustable |
| Return Speed | 20 mm / s, fixed |
| Vibration Frequency | 0-20 Hz, adjustable |
| Vibration Amplitude | 0.3 mm, fixed |
| Z-axis Vibration | ~0 μ m |
| Blade | ceramic blades |
| Cutting Angle | 20 degrees, fixed (standard) |
| Thickness Adjustment | digital adjustment, automatic |
| Minimum Slice Thickness | 10 μ m |
| Maximum Adjustable Thickness | 250 μ m |
| Slice Thickness Resolution | 1 μ m |
| Maximum Tissue Diameter | 22 mm |
| Maximum Tissue Length | 25 mm |
| Cutting Mode | single/continue selectable |
| Bath | 210 x 80 x 35 mm |
| Power Source | DC 13-15 V |
| Power Consumption | 4 W |
| Dimension (L x W x H) | 300 x 280 x 280 mm |
| Weight | 9kg |

Precisionary Instruments LLC
90 Canal Street, Suite 400
Boston, MA 02114

Tel: 617-682-0586
Fax: 1(866)424-2217
Email: info@precisionary.com