

VF-700-0Z

Technical specs



PRECISIONARY



Advance Speed	0-20 mm / s, adjustable
Return Speed	20 mm / s, fixed
Vibration Frequency	0-20 Hz, adjustable
Vibration Amplitude	0.3 mm, fixed
Z-axis Vibration	~0 μ m
Blade	ceramic blades
Cutting Angle	20 degrees, fixed (standard)
Thickness Adjustment	digital adjustment, automatic
Minimum Slice Thickness	10 μ m
Maximum Adjustable Thickness	250 μ m
Slice Thickness Resolution	1 μ m
Maximum Tissue Diameter	22 mm
Maximum Tissue Length	25 mm
Cutting Mode	single/continue selectable
Bath	210 x 80 x 35 mm
Power Source	DC 13-15 V
Power Consumption	4 W
Dimension (L x W x H)	300 x 280 x 280 mm
Weight	9kg

Precisionary Instruments LLC
207 Union Street
Natick, MA 01760

Tel: 617-682-0586
Fax: 1(866)424-2217
Email: info@precisionary.com