



# COMPLETE TISSUE EMBEDDING CENTER TEC-200

## User Manual

Includes:  
Paraffin Dispenser  
Cold Plate

[www.precisionary.com](http://www.precisionary.com)



Abigail Chu, MD-PhD  
Chief Scientific Officer (CSO)  
Precisionary Instruments  
617-682-0586

## Technical Specifications

<b><i>Embedding Center</i></b>	Paraffin chamber capacity	5.5. L
	Maximum temperature	75°C
	Temperature control precision	±1%
	Paraffin flow control	Touch screen and optional foot pedal
	Working voltage	AC 220V+10% 50Hz (standard model); AC 110V+10% 60Hz
	Power	1000W
	Dimensions	620x560x410mm (W x D x H)
	Net weight	39kg
	Dimensions of paraffin chamber	265x145x145mm (W x D x H)
<b><i>Cooling Plate</i></b>	Minimum temperature	-24°C
	Suggested working temperature	-10°C
	Working voltage	AC 220V+10% 50Hz (standard model); AC 110V+10% 60Hz
	Dimensions	620x380x410mm (W x D x H)
	Net weight	25kg
	Dimensions of cryo module	345x330mm