



Precisionary CF-6100 Cryostat User Manual

www.precisionary.com



Abigail Chu, MD-PhD
Chief Scientific Officer (CSO)
Precisionary Instruments
617-682-0586

Technical Specifications

<i>Parameter</i>	<i>Specification</i>
Nominal supply voltage	220V~240V/110V
Nominal frequency	50Hz/60Hz
Main fuse	4A
Nominal power	1200VA
Operating temperature range	+5°C to +30°C
Transport and storage temperature	+5°C to +55°C
Relative humidity	Maximum 60%, non-condensing
A weighted noise level	≤ 70 dB (A)

Size and Weight

"Size, without packaging, excluding hand wheel (width × depth × height)"	650 × 700 × 1280mm, ±5%
"Size, without packaging, including hand wheel (width × depth × height)"	780 × 700 × 1280 mm, ±5%
Weight, without packaging	140kg
Weight, with packaging	170kg